


## How to Use

### Home-care Beauty Device Liftlogy YJFD1N

#### Names of Components

##### Main unit



##### Accessories

- USB charging cable (USB Type-C)
- Skincare product
- User's Manual

#### Modes

| MODE indicator | Photo      | Areas of use                     | WARMTH Level | EMS Level | Recommended duration of use for each treatment | Recommended frequency of use |
|----------------|------------|----------------------------------|--------------|-----------|--|------------------------------|
| FACE           | Blue green | Face and neck                    | 1-3          | 1-6       | 2 minutes<br>(Max 10 minutes)                  | Every day                    |
| BODY           | Red        | Decollete, upper arms, and flank |              |           |  |                              |

##### Recommended use

The device can be used every day. However, as a general guide, please use the device once a day, up to 2 minutes per treatment (max 10 minutes on each area) for a total of up to 20 minutes a day.

Have an interval of at least 24 hours when using again the following day and beyond.

\*This is only a guideline. Please observe your skin condition and use the device accordingly. Excessive use may cause skin damage.

##### Precautions

Read the items below carefully before use and use only as directed.

- Remove all metal accessories (watches, necklaces, rings, earrings, etc.) when using the device.

Failure to comply may cause an accident or an injury.

- Be sure to remove makeup before using the device.

If you use the device over makeup, the foundation may cause dark staining of the skin.

- Make sure to perform a pretest before using the device for the first time.

- Moisten your skin with the included skincare product or lotion, gel, or water before applying directly to your skin.

If using your regular skincare product, please use a thick lotion or serum.

## How to use

### 1 Turn on the device

Press and hold the POWER/WARMTH LEVEL selection button (approx. 1 second).

### 3 Select the WARMTH level

Short press the POWER/WARMTH level selection button to select the warmth level from levels 1 to 3.

### 5 Place the device on your skin

Gently place the main treatment head on your skin.

### 2 Select the mode

Short press the MODE selection button to select the mode.

### 4 Select EMS level

Short press the EMS level selection button to select the EMS level from levels 1 to 6.

### 6 Turn off the device

Press and hold the POWER/WARMTH LEVEL selection button (approx. 2 seconds).

## How to move the device over each area

Moving the main treatment head as follows is effective.

### FACE Mode

Recommended usage: 1 set (2 minutes) \*Max 5 sets (10 minutes)

Please use by pressing in the order of ① to ⑥ as shown in the figure.

- Use for 1 minute for half-face, 2 minutes for full face.  
When the half-face alarm sounds after 1 minute, switch to the other side.  
If 1 set (2 minutes) is not enough, use again in the order of ① to ⑥.

- ① Neck: hold for 10 seconds
- ② Lower jawline (below face line): hold for 10 seconds
- ③ Cheeks\*: hold for 20 seconds  
\*The recommended position is on the line connecting the corner of the mouth and the ear.
- ④ Temple: hold for 5 seconds
- ⑤ Edge of forehead: hold for 5 seconds
- ⑥ Cheeks: gently slide the device on skin surface\* for 10 seconds  
\*Hold the device perpendicular to the mouth and use with the image of moving diagonally upward.




- The device will beep after 1 minute of selecting a mode (timer function).  
After another 1 minute (2 minutes in total), the device will beep three times and enter standby status (standby function).

### BODY Mode

Recommended usage: 1 set (2 minutes) \*Max 5 sets (10 minutes)


#### Decollete

Move from inside to outside along the collarbone.



#### Upper arms

Move from the elbow towards the armpit.



## ATTENTION

- The EMS function may stimulate the nerves in your teeth. Do not press the main treatment head strongly against your skin but instead lower the level or move the position. Be especially careful if you are undergoing dental treatment or have dental fillings.
- When using the device on your neck, be careful not to press too hard, and stop use immediately if you feel unwell during use. Be careful as you may feel stronger EMS stimulation on some occasions.